


Improve business efficiency • Improve employee workplace safety • Improve patient care & management



#153055
w/Lavender
#153046

#153055B w/Norway #153035

Beamic™ Screen Frame
#153055 - #153055B

The 5 sections of the screen are attached to each other by a riveted hinge which allows the frame to smoothly fold and unfold as needed. The S-shaped base of each screen section on the frame also increases stability for use in busy environments. Privacy screen panels, either in economical vinyl or attractive designer fabrics are needed to complete the screen. One panel is required for each section and is secured by top and bottom snap in rods. Panels are discounted when purchased with the screen.

Features:

- Available in 2, 3, 4 or 5 section models
- Five sections measures 68”H x 91”W
- Requires one vinyl or fabric panel per section
- Anodized aluminum construction provides a light weight and durable solution eliminating the possibility of paint chipping and rust
- 5 year warranty on Frames
- 90 Days on Screen Panels




White
#153010



Green
#153015

Screen Frames and Panels are Sold separately




Apricot
#153020




Oatmeal
#153028



Marine
#153029



Norway
#153035




Sage
#153048



Rouge
#153040




Autumn
#153042



Powder
#153044



Frost
#153045



Lavender
#153046