

Improve business efficiency • Improve employee workplace safety • Improve patient care & management



#153961 w/Rouge - #153040

3-Section King Size Screen Frame with casters #153961

3 Section Mobile King Screens are made of a durable, lightweight and easy to clean anodized aluminum. These screens are secured to a heavy duty cross linked U shaped base on four 2 inch casters for easy mobility and increased stability. The 3 sections of the screen are attached to each other by a riveted hinge which allows the frame to smoothly fold and unfold as needed. Privacy screen panels, either in economical vinyl or attractive designer fabrics are needed to complete the screen. Two panels are required for each section and are secured by top and bottom snap in rods. Panels are discounted when purchased with the screen. King screens are the largest of the Omnimed privacy screens and are ideal for areas where maximum privacy is a must.

Features:

- Frame utilizes casters for quiet and effortless mobility
- Three sections measures 73 $\frac{3}{4}$ "H x 96"W
- Requires two vinyl or fabric panels per section
- 5 year warranty on Frames
- 90 Days on Screen Panels




White
#153010



Green
#153015

Screen Frames and Panels are Sold separately




Apricot
#153020




Oatmeal
#153028



Marine
#153029



Norway
#153035




Sage
#153048



Rouge
#153040




Autumn
#153042



Powder
#153044



Frost
#153045



Lavender
#153046