



## **Narcotic Cabinet with Combination Lock - #181680**

Constructed of durable, heavy duty stainless steel and meet schedule 2 drug storage requirements. They have strong piano hinged doors and a two door, dual locking system to ensure content security. The outside door has a triple bolt push button combo lock which opens to a second door secured by a second combo lock with a heavy gauge locking cam. The ambi top of the cabinet lifts so that the unit can be flipped to accommodate left or right outer door access. Its four aluminum shelves are adjustable and increase storage capacity. Pre-drilled holes that are reinforced with stainless steel mounting discs ensure stable wall attachment. These units are warranted for life & include 1 year on the locks.

- Stainless steel double door cabinet with 4 shelves
- DoubleDoor Triple bolt push button combo lock
- Ambi top for left/right outer door access
- Pre-drilled holes reinforced with stainless steel mounting discs
- Measures 24”H by 16”W by 8”D
- Unit Has lifetime Warranty (1 Year On Locks)