

OMNIMED Specimen Dropbox Cabinet

Improve business efficiency • Improve employee workplace safety • Improve patient care & management



Specimen Dropbox Cabinet #181790 & 181790T

Constructed of durable, heavy duty stainless steel. Provides a convenient way to secure specimens until a clinician can retrieve. This unit is wall mounted so it will not affect the fire rating like a traditional pass-thru cabinet. These units are warranted for life & include 1 year on the locks.

- 3 Independent specimen chutes
- Holds up to 9 specimen cups
- 181790 has key lock security - Multiple units keyed the same unless specified
- 181790T has thumb -latch for easy access by patient or clinician
- Unit Has lifetime Warranty (1 Year On Locks)
- Overall Dimensions: 10"H X 13.5"W X 3.5" D
Chute Dimensions: 9.5"H X 2.5"W X 2.5"D

OMNIMED

800.257.2326 • www.omnimedstore.com