

National 205

Compact yet Powerful

National's 205 is an ideal microscope for teaching microscopy to primary levels or for the hobbyist. Model 205 is a monocular microscope that supports an 18mm field of view, similar to larger and more expensive models.

FEATURES:

Easy grip design to enhance easy handling within classrooms

Low position focusing controls

State of the Art LED illumination

Patented 1-touch stage clip enables easy change of slides

Inward facing nosepiece provides maximum operating space

- 185 mm Focal Length
- 45° inclined monocular
- WF10X/18mm eyepieces
- 4x,10x,40x objectives
- Separate coarse and fine focus
- 70mW/3.5v LED illumination w/intensity control



■ National 205