

Computer Monitor Stand #350300**Construction:** Plated steel and aluminum

Omnimed computer monitor stands have a stainless steel Laptop tray secured to a heavy gauge plated steel column. The unit's easy grab wheel ring handle comes with one Poly tray and can accommodate up to four trays within the ring. The sturdy 5 leg base and 3 inch casters allows for easy transport. These stands can be outfitted with a variety of accessories including transport baskets of various sizes, a printer tray and additional Poly trays for the wheel ring. The work surface measures 16 9/16W x 13 1/6L and stands 42 inches from the floor. The base has a diameter of 25 inches. These units include a 5 year warranty.



Computer Monitor Stand
#350300
Shown with optional:
Wheel Ring Handle - #350051
Poly Tray - #350055



Wheel Ring Handle
#350051



Poly Tray
#350055



Wire Basket
6 1/4H x 4 1/2W x 4 1/4D
#350005



Transport Basket
10H x 13W x 4D
#350002