

*Improve business efficiency • Improve employee workplace safety • Improve patient care & management*

---

## **Suture Transport Cart #350309**

The Suture Transport Cart is designed to mobilize a large assortment of sutures. The small platform is perfect for doctor's offices, urgent care facilities, emergency rooms and hospital floors. The stand height can easily be adjusted from 31" to 41" so that the sutures can be accessed from a seated or standing position.



- 18 suture bins per cart
- Made from impact resistant acrylic
- Height adjustment of tray is 31" to 41"
- Dimensions of tray with handles – 21.5" W x 14.5" D
- Stainless steel tray to support suture bins
- Aluminum grab and go handles
- 22" diameter star base
- 1-Year Warranty

## **Triple Suture Rack #304102**

The Triple Suture Rack is made of special impact resistant acrylic that is stronger and longer lasting than traditional acrylic. It is easy to keep clean with simple soap and water.

- 9 bins per rack
- Open in front for easy restocking
- Mounts to walls, countertops or transport stands
- Dimensions: 15.75" W x 5.5" D 18.25" H

