

Improve business efficiency • Improve employee workplace safety • Improve patient care & management

All computing Work Stations keep patient information concealed and computers secure.
Our cost-effective innovative designs are user friendly, durable and ergonomically efficient.



Laptop / Vital Signs Cart
#350750

Laptop / Vital Signs Cart #350750

Features:

Color: White / Custom colors also available

Construction: Powder coated heavy gauge steel

- Low center of gravity base with tough powder-coat finish
- Easy rolling 4" rubber casters
- Telescoping column with pneumatic assist for easy lifting
- 46" fully extended and 36" compressed
- Left to right sliding mouse tray
- Easy access ring style handle for easy grab and go use
- Fully adjustable 19" x 7" keyboard tray
- One poly tray in handle for additional storage
- Removable basket - 8"H x 10"W x 4"D
- Laptop security head
- Scanner and mouse holder
- Vital Signs monitor holder
- 10' hospital grade retractable cord reel
- 5-year warranty