



## Power Lifters

### **Blood Filtration Power Lifter #741318**

Single Column System

**Features:**

Safe air-operated system featuring an internal, rechargeable 12V battery

Easy to operate:

- Foot pedal operates pneumatic air-glide system
- Column extends from 5' 7" to 7' 5" high
- Accommodates use of various leuko-depletion filters
- 18 powder coated steel hooks
- 6 removable blood collection sections. Each section has 3 bins for a total of 18 bins
- Rechargeable 12V battery & charger
- 6 non-conductive heavy-duty casters
- Optional - Low Voltage Light Kit #741316
- Warranty- 5 Years

