

# Palm Abbe™

## Digital Refractometer

### BRIX REFRACTOMETER

The MISCO PA201 digital Brix refractometer is a valuable instrument for determining the sugar content in fruits, vegetables, juices, beverages, and in wine and beer making. It is also used extensively as a base scale for measuring non-sugar fluids such as metalworking fluids, lubricants and industrial chemicals. Range 0 to 56 Brix, Resolution 0.1%, Precision +/-0.1%.

Model: PA201 Part No.: PA201

Scale #	Unit of Measure	Range	Resolution	Precision (+/-)
004	Brix	0 to 56	0.1	0.1



Image: PA202 shown with optional rubber armor jacket.



Shown with optional rubber armor jacket.

	MISCO Palm Abbe	Competitor
Prism Material	SAPPHIRE (Comparable to Diamond)	GLASS
Detector Array Elements	1,024 Elements	128 Elements
Detector Resolution	3,256 ppi	400 ppi
Stainless Steel Well	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Number of Scales	Up to Five	One
Display	24 Alphanumeric Characters x2 Lines	Four Digits
Display Backlight	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Evaporation Cover	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Optional Rubber Armor	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Made In	U.S.A	Asia

Made in the USA! - Copyright © 2013 MISCO



www.misco.com

3401 Virginia Road  
Cleveland, Ohio 44122 USA  
TEL. 216-831-1000  
Toll Free 800-358-1100