

Palm Abbe™

Digital Refractometer

BRIX & REFRACTIVE INDEX REFRACTOMETER

Dual-scale digital refractometer with scales for measuring both Brix (to 85) and refractive index (to 1.5000). Brix is a measure of sugar concentration and is used extensively in the agriculture and food and beverage industries. Refractive index is a fundamental physical property of a solution and is related to the speed light travels through it. Both the Brix and refractive index scales can be used to measure a nearly unlimited number of fluids.

Model: PA202 Part No.: PA202

Scale #	Unit of Measure	Range	Resolution	Precision (+/-)
001	Refractive Index (nD)	1.3330 to 1.5000	0.0001	0.0001
003	Brix	0 to 85	0.1	0.1



Image: PA202 shown with optional rubber armor jacket.



Shown with optional rubber armor jacket.

	MISCO Palm Abbe	Competitor
Prism Material	SAPPHIRE (Comparable to Diamond)	GLASS
Detector Array Elements	1,024 Elements	128 Elements
Detector Resolution	3,256 ppi	400 ppi
Stainless Steel Well	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Number of Scales	Up to Five	One
Display	24 Alphanumeric Characters x2 Lines	Four Digits
Display Backlight	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Evaporation Cover	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Optional Rubber Armor	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Made In	U.S.A	Asia

Made in the USA! - Copyright © 2013 MISCO



www.misco.com

3401 Virginia Road
Cleveland, Ohio 44122 USA
TEL. 216-831-1000
Toll Free 800-358-1100