MD300 C2 Pulse Oximeter

SENSING

ANALYSIS

DELIVERY



MD300 C2 PULSE OXIMETER OVERVIEW

The ChoiceMMed MD300 C2 is a great addition to the line of pulse oximeters offerd by Maxtec. Its compact design allows it to be used for adult applications in many different areas such as hospital, homecare, and EMS. This product is an excellent source for the accurate pulse oximeter that you've been looking for.

PART NUMBERS

MD300 C2 Pulse Oximeter.....R204P23

Specifications

Adult-Pediatrics Patient Range

SpO₂

70% - 100% Measurement Range

Resolution 1%

Measurement Accuracy 70% - 100%: ± 3%; 70%: No Definition

Pulse Rate

Measurement Range 30 - 235 bpm

Resolution 1 bpm

T 30 - 99: ± 2%; 100 - 235: ± 2% Measurement Accuracy

Display

Dual Color OLED Display Туре

SpO₂, PR, Pulse Bar, Plethysmogram Parameters

6 Display Modes Mode Adjustable 10 Levels **Brightness**

Alarm

Alarm **Battery Low Indicator**

Mode Visual

Mechanical

58mm (L) x 32mm (W) x 37mm (H) Dimensions

Weight 33g (without batteries)

Battery

Battery: Two AAA 1.5V, 800mAh alkaline batteries (approx. 30 hours)

Environmental

5°C - 40°C Operation Temperature Storage Temperature -20°C - 55°C

Operation Humidity <80%, no condensation Storage Humidity <93%, no condensation

This product complies with IEC60601-1-2 for electromagnetic compatibility Class $\ensuremath{\mathsf{B}}$ Electromagnetic Compatibility

Interference Resistance Capacity Device works normally when mixed noise produced by BIO-TEK INDEX

Against Ambient Light Pulse Oximeter tester

