MD300 C63 Pulse Oximeter



MD300 C63 PULSE OXIMETER OVERVIEW

The **MD300 C63** spot-check pulse oximeter provides exceptional quality and performance in a small, robust design. This model includes new advanced features like a shockprook design, finger skid resistance and automatic compensation for infrared signal interference. This unit is suitable for use in hospitals, home care and even during exercise. PART NUMBERS MD300 C63 Pulse Oximeter......R204P27

Specifications

Patient Range

SpO₂

Measurement Range Resolution Measurement Accuracy

Pulse Rate

Measurement Range Resolution Measurement Accuracy

Display

Type Parameters Mode Brightness

Alarm

Alarm Mode

Mechanical

Dimensions Weight

Battery

Battery:

Environmental

Operation Temperature Storage Temperature Operation Humidity Storage Humidity Electromagnetic Compatibility

Interference Resistance Capacity Against Ambient Light Adult-Pediatrics

70% - 100% 1% 70% - 100%: ± 3%; 70%: No Definition

30 - 235 bpm 1 bpm T 30 - 99: ± 2%; 100 - 235: ± 2%

Dual Color OLED Display SpO₂, PR, Pulse Bar, Plethysmogram 6 Display Modes Adjustable 10 Levels

Battery Low Indicator Visual

58mm (L) x 32mm (W) x 37mm (H) 33g (without batteries)

Two AAA 1.5V, 800mAh alkaline batteries (approx. 30 hours)

5°C - 40°C -20°C - 55°C <80%, no condensation <93%, no condensation This product complies with IEC60601-1-2 for electromagnetic compatibility Class B Device works normally when mixed noise produced by BIO-TEK INDEX Pulse Oximeter tester



ML# 224 REV. E