


# SF, Stainless Steel Filter Holder



## Features

- ▶ Spin-lock connection (M577346) facilitates fast and stable installation without clamp.
- ▶ Laser-marked graduation helps accurately sample addition.
- ▶ Membrane guide on the support base helps position membrane easily.
- ▶ Forceps access points on support base help retrieve membrane without damage.
- ▶ SS316 funnel is shatterproof, durable and can be fast sterilized by flame.

## Specifications

Material	100 ml	300 ml	500 ml
Funnel	SS316	SS316	SS316
Support base	SS316	SS316	SS316
Membrane support	SS316	SS316	SS316
Funnel lid	SS316	SS316	SS316

Specifications	100 ml	300 ml	500 ml
Funnel capacity	100 ml	300 ml	500 ml
Filter diameter	47 mm	47 mm	47 mm
Effective filtration area	11.6 cm <sup>2</sup>	11.6 cm <sup>2</sup>	11.6 cm <sup>2</sup>
Graduated increment	50 / 100 ml	100 / 250 / 300 ml	250 / 500 ml

## Ordering Information

<b>SF 1, Stainless Steel Filter Holder (100 ml)</b>		<b>180100-10</b>
1	Stainless steel funnel, 100 ml	180100-01
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16
<b>SF 3, Stainless Steel Filter Holder (300 ml)</b>		<b>180100-30</b>
1	Stainless steel funnel, 300 ml	180100-03
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16
<b>SF 5, Stainless Steel Filter Holder (500 ml)</b>		<b>180100-50</b>
1	Stainless steel funnel, 500 ml	180100-05
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16
<b>SF, Stainless Steel Lid</b>		
	100 ml Stainless steel lid	180100-11
	300 ml / 500 ml Stainless steel lid	180100-13

