



Hot/Cold Therapy & Clinic Supplies



## Introduction

Chattanooga® is a leading manufacturer of rehabilitation equipment for treating musculoskeletal, neurological, and soft tissue disorders. For over seven decades, Chattanooga products have helped contribute to positive treatment outcomes in hospitals, clinics, and home settings worldwide. We strive to lead the rehabilitation industry by example, and through continuous innovation, provide real world solutions for clinicians and their patients.

As one of the leading brands in rehabilitation, Chattanooga is sold in over 80 countries worldwide. Chattanooga dealers, distributors and customers are supported by a seasoned team of personnel who oversee product development, product training and education, order placement, product delivery, after sales service, and provide dedicated sales and marketing support.

Chattanooga is a brand of DJO Global. DJO Global provides solutions for musculoskeletal health, vascular health, and pain management. Our products help prevent injuries, rehabilitate after surgery, injury or from degenerative disease, helping to enable patients to regain or maintain natural motion.

Please read complete instructions for use for all products prior to use.

Visit [www.DJGlobal.com/chattanooga](http://www.DJGlobal.com/chattanooga)



AIRCAST

 chattanooga

CMFV

Compex

djosurgical

DONJOY

Dr.  
Comfort

PROCARE

## Contents

Hydrocollator® Mobile Heating Units	2	Stabilizer™	10
Hydrocollator® Stationary Heating Units	3	Sully® Shoulder Support	11
Hydrocollator® HotPacs™	4	Gels & Lotions	12
Hydrocollator® Terry Covers	5	Electrodes & Accessories	13
ColPac® Chilling Unit	6		
ColPac®, Boo-Boo Pac™, Flexipac™, Sport-Pac™	7		
Theratherm Digital Moist Heat	8		
Carts & Stools	9		



## Hot/Cold Therapy and Clinic Supplies

A leader in therapeutic products, Chattanooga offers a diverse selection of high quality hot/cold therapy and clinical rehabilitation supplies. Chattanooga is the proud inventor of the original Hydrocollator HotPacs and heating unit, a product that was first to market more than 70 years ago, and is still widely used by physical therapists and other medical professionals when using heat therapy. The HotPac generates deep moist heat which helps relieve chronic pain, acute muscle strains and general musculoskeletal discomfort. Chattanooga also raised the bar in cold therapy with the original ColPac, which comes in a variety of shapes and sizes. Chilling units are also available for convenient storage and cooling. For athletic training rooms and physical therapy clinics; carts, stools, lotion, and electrodes can all be ordered. Combining industry leading hot and cold therapy products with clinic supplies, Chattanooga can meet many clinics needs. Please read all instructions for use and any additional materials prior to use.

# Hydrocollator® Mobile Heating Units

## Ordering Information

Item#	Description
-------	-------------

### M-2 Mobile Hydrocollator

2402	M2 Mobile - Includes 1 Standard (1006) HotPacs
402-2	M2 Mobile - Includes 3 Standard (1006), 3 Oversize (1004), and 3 Cervical (1002) HotPacs
2402-3	M2 Mobile - Includes 4 Standard (1006), 2 Oversize (1004), and 4 Cervical (1002) HotPacs
2402-4	M2 Mobile - Includes 4 Standard (1006), and 4 Oversize (1004)

### M-4 Mobile Hydrocollator

2502	M4 Mobile - Includes 24 Standard (1006) HotPacs
2502-2	M4 Mobile - Includes 8 Oversize (1004), and 8 Cervical (1002) HotPacs
2502-3	M4 Mobile - Includes 6 Standard (1006), 6 Oversize (1004), and 6 Cervical (1002) HotPacs

### SS-2 Mobile Hydrocollator

2302	SS2 Mobile - Includes 8 Standard (1006) HotPacs
2302-2	SS2 Mobile - Includes 4 Standard (1006), and 4 Cervical (1002), HotPacs
2302-3	SS2 Mobile - Includes 2 Standard (1006), 2 Oversize (1004), and 2 Cervical (1002) HotPacs

## Technical Specifications

### M-2 Mobile Hydrocollator

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	48 lb (22 kg) /Shipping - 93 lb (42 kg)
Dimensions*:	27" x 16" x 33" (69 x 41 x 84 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 6011

### M-4 Mobile Hydrocollator

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1500 W
Weight:	132 lb (60 kg) /Shipping - 160 lb (72 kg)
Dimensions*:	35" x 20" x 33" (89 x 51 x 84 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 6011

### SS-2 Mobile Hydrocollator

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	66 lb (30 kg) /Shipping - 80 lb (36 kg)
Dimensions*:	21" x 16" x 33" (53 x 41 x 84 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 6011

\* Dimensions are listed as L x W x H



Our benchmark product is durable and easy to maintain, these high-quality stainless steel units give you a constant supply of temperature-consistent HotPacs.

## Features

### Hydrocollator Mobile Heating Units

- High-quality stainless steel
- Thermostatically controlled temperature with safety feature
- 3" (8 cm) swivel-type rubber casters for silent, friction-free movement
- Easy maintenance. Simple to fill and drain. No plumbing required
- One-year warranty

\* Approximate with HotPacs

	M-2	M-4	SS-2	E-2	E-1/E-1 MHT
Tank Capacity	14 gal (52 L)	27 gal (102 L)	10 gal (37 L)	8 gal (37 L)	3 gal (37 L)
Temperature Range	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)
Thermal Cut-Out Temp.	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)
Temperature Accuracy	+/- 10%	+/- 10%	+/- 10%	+/- 10%	+/- 10%
Heat up Time to 160° F (to 70°C)	6 Hours	8 Hours	4 Hours	4 Hours	2 Hours
Cool Down Time from 160° F (from 70°C)	3 Hours	4 Hours	4 Hours	2 Hours	1 Hour
Fiberglass Insulation	Yes	Yes	Yes	No	No

Be sure to read the IFU for all warnings and precautions

## Hydrocollator® Stationary Heating Units



## Hydrocollator Accessories



### Wall Mounted Towel Rack

Six hooks extend 13" (33 cm) from wall. Measures 2" x 33" (5 cm x 84 cm).



### Mini Non-Electric Lotion Warmer

Hangs on the front of the Hydrocollator E-1 or E-2 heating unit and uses the unit's warmth to keep lotion or gel at a comfortable temperature. Bottles are not included.



### TM-1 Electric Lotion Warmer

Easily attaches to any treatment table or wall to keep massage lotion or ultrasound gel at a comfortable temperature and within reach. Bottles are not included.

Be sure to read the IFU for all warnings and precautions

### Ordering Information

Item#	Description
E-1 Stationary Hydrocollator	
2102	E1 Stationary - Includes 4 Standard (1006) HotPacs
2102-2	E1 Stationary - Includes 3 Standard (1006) HotPacs and 1 Cervical (1002) HotPac
2102-2	E1 Stationary - Includes 2 Standard (1006) HotPacs and 2 Cervical (1002) HotPacs

### E-2 Stationary Hydrocollator

2802	E1 Stationary - Includes 2 Oversize (1004), 3 Standard (1006) HotPacs and 1 Cervical (1002) HotPac
------	--

### Hydrocollator Accessories

4016	Wall Mounted Towel Rack
22100	Mini Non-Electric Lotion Warmer
4012	TM-1 Electric Warmer
4228	Dial Thermometer (not shown)

### Technical Specifications

E-1 Stationary Hydrocollator	
Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	15 lb (7 kg) /Shipping - 21 lb (10 kg)
Dimensions*:	13" x 8" x 16" (33 x 20 x 41 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

E-2 Stationary Hydrocollator	
Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	29 lb (13 kg) /Shipping - 46 lb (21 kg)
Dimensions*:	15" x 13" x 20" (38 x 33 x 51 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

\* Dimensions are listed as L x W x H

**Ordering Information**

Item#	Description
Hydrocollator HotPac - Moist Heat Hot Pack	
A 1008	10" x 18" (25 cm x 46 cm)
B 1010	10" x 24" (25 cm x 61 cm)
C 1004	Oversize 15" x 24" (38 cm x 61 cm)
D 1002	Neck Contour 24" long (61 cm)
E 1014	Half Size 5" x 12" (13 cm x 30 cm)
F 1006	Standard 10" x 12" (25 cm x 30 cm)
G 1012	Knee-Shoulder 10" x 20" (25 cm x 50 cm)

# Hydrocollator® Hotpacs™



A leading choice of clinicians for over 70 years

**Features**

**The Original Hydrocollator® Moist Heat HotPac**

Chattanooga is the inventor of the original Hydrocollator HotPac and proud to still be producing this product today.

**Hydrocollator HotPacs**

- Seven HotPac sizes designed to fit every body contour
- Provides up to 30 minutes of deep, soothing, moist heat
- Brown tabs on most HotPacs for easy rotation in busy clinics
- One HotPac provides hundreds of cost effective treatments
- One-year warranty



# Hydrocollator® Terry Covers



## Features

### All-Terry Covers

- For high-volume institutional use or therapy at home
- Hook and loop closures to prevent unwrapping or slipping
- Washable and reusable
- One-year warranty

### Foam-Filled Terry Covers

- For home or light clinic use
- Polyurethane foam layer between two layers of terry
- Available with or without pockets
- Washable and reusable
- One-year warranty



## Ordering Information

Item#	Description
<b>All-Terry Covers</b>	
1118	Standard 19" x 27" (48 cm x 69 cm)
1126	Oversize 24" x 36" (61 cm x 91 cm)
1104	Neck Contour 25" x 18" (64 cm x 46 cm)
<b>Foam-Filled Terry Covers</b>	
1102	Standard 19" x 27" (48 cm x 69 cm)
1124	Oversize 24" x 36" (61 cm x 91 cm)
1120	Neck Contour 25" x 16" (64 cm x 41 cm)
<b>Foam- Pocket Terry Covers</b>	
1108	Standard 17" x 27" (43 cm x 69 cm)
1110	Oversize 24" x 36" (61 cm x 91 cm)
<b>Hydrocollator HotPac Set in Retail Packaging</b>	
1064	Standard Set in Retail Package includes one Standard HotPac and one Foam Filled Terry Cover
1062	Neck Set in Retail Package includes one Neck HotPac and one Foam Filled Terry Cover

# ColPac® Chilling Unit



## Ordering Information

Item#	Description
-------	-------------

### ColPac Chilling Unit Model C-5

3102	C-5 ColPac Chilling Unit Includes 6 Standard (1500) and 6 Half Size (1506) ColPacs
------	--

### Standard Accessories

4010	Side Table Rack See page 4 for Description and image
------	--

4230	Extra Shelf for Side Table Rack (not shown)
------	---

### ColPac Freezer

90910	ColPac Freezer 21 1/4" x 26" x 33 1/2" (53 cm x 66 cm x 86 cm) Includes 12 Standard ColPacs (1500) Shipping weight 150 lbs (68 kg)
-------	--

## Technical Specifications

### ColPac Chilling Unit Model C-5

Mains Power:	115 VAC- 60 Hz
--------------	----------------

Current:	4A
----------	----

Weight:	165 lb (74 kg) Shipping Weight - 190 lb (86 kg)
---------	---

Dimensions*:	27" x 16" x 33" (69 x 41 x 84 cm)
--------------	-----------------------------------

Electrical Safety Class:	Type B
--------------------------	--------

Safety Tests:	UL 544 (U.S. Only)
---------------	--------------------

\* Dimensions are listed as L x W x H



3102



90910

Large capacity and fast chilling time means your patients will never have to wait for a ColPac. The ColPac Chilling Units help you keep pace throughout a busy day. No plumbing required, simply plug into a standard electrical outlet.

### Features

#### ColPac Chilling Unit

- Closed-cell foam insulation for efficiency
- Stainless steel, welded construction for long life
- Heavy-duty compressor for cooling
- Specially designed interior chilling coils cool packs fast
- 3-inch (8 cm) swivel-type rubber casters for silent, friction-free movement of mobile units
- Drain valve for easy cleaning and defrosting
- One-year warranty

#### ColPac Freezer

- Five cubic feet of spacious capacity
- Four durable stainless steel shelves
- Comes complete with trays and cup holders
- Shipped with 12 standard size ColPacs
- MetLabs Certified in the US and Canada
- One-year warranty





# ColPac®, Boo-Boo Pac™, Flexipac™, Sport-Pac™

Seven different sizes of blue vinyl and four sizes of black polyurethane ensure cold therapy for any body area or shape.



## Features

### The Original ColPac

- Designed for heavy clinical use.
- Provides pain management
- Filled with a silica gel that remains pliable throughout treatment
- Remains cold for up to 30 minutes
- Not made with natural rubber latex
- One-year warranty

### boo-boo pac

- Fun, bear-shaped design especially made for kids
- Ideal for pediatric applications
- Filled with silica gel
- Remains cold for up to 30 minutes
- Shipped in a display box containing 12 packs

### Flexipac

- Freezable to provide penetrating cold therapy
- Microwavable for soothing heat
- Flexible gel conforms to any body contour at any temperature

### Sport-Pac

- Soccer ball design is ideal for youth sport applications
- Easily transferred from the freezer to the team cooler for sideline use
- Remains cold for up to 30 minutes
- Shipped in a display box containing 10 packs

*These products can be cold enough to cause serious injury. Read and understand all warnings and instructions for use before using the product.*

## Ordering Information

Item#	Description
<b>Blue Vinyl ColPac</b>	
A 1512	Over Size 11" x 21" (28 cm x 53 cm)
B 1500	Standard Size 11" x 14" (28 cm x 36 cm)
C 1506	Half Size 7½" x 11" (19 cm x 28 cm)
D 1508	Neck Contour 23" long (58 cm)
E 1502	3" x 11" (8 cm x 28 cm)
F 1504	Quarter Size 5½" x 7½" (14 cm x 19 cm)
G 1510	Eye (Specially shaped to cover closed eyelids)
<b>Black Polyurethane ColPac</b>	
H 1554	Neck Contour 21" long (53 cm)
I 1556	Over size 12½" x 18½" (32 cm x 47 cm)
J 1552	Standard Size 10" x 13½" (25 cm x 34 cm)
K 1562	Half Size 6½" x 11" (17 cm x 28 cm)
<b>L boo-boo pac</b>	
1534	Royal Blue boo-boo pac (12 per case) Shipped in a display box containing 12 packs. NOTE: Must be purchased in quantities of 12 each
<b>M FlexiPac Hot &amp; Cold Compresses</b>	
4026	5" x 6" (13 cm x 15 cm) 48/case
4020	5" x 10" (13 cm x 25 cm) 24/case
4029	8" x 14" (20 cm x 35 cm) 12/case
4032	5" x 10" (13 cm x 25 cm) 25/case Chiropractically Yours (not shown)
<b>N Sport-Pac</b>	
1517	Sport-Pac (10 per case) Shipped in a display box containing 10 packs. NOTE: Must be purchased in quantities of 10 each

## Ordering Information

## Item# Description

## Theratherm Digital Moist Heat Pack

1030	Small 7" x 15" (18 cm x 38 cm)
1031	Medium 14" x 14" (36 cm x 36 cm)
1032	Standard 14" x 27" (36 cm x 68 cm)
1033	Shoulder/Neck 23" x 20" (58 cm x 50 cm)

## Theratherm® Digital Moist Heat Pack



Theratherm Digital Moist Heat Pack provides moist heat therapy for pain management. Its flannel cover draws moisture from the air and as the pack's temperature rises, the moisture is released onto the skin, speeding relief to your patients. Its unique, easy-to-use digital control is convenient.

## Features

## Theratherm Digital Moist Heat Pack

- Convenient preset time and temperature
- Digital hand control for operator to program and monitor temperature and/or treatment time
- Lockout mode protects against inconsistent setting changes throughout treatment sessions
- Moist heat therapy temperatures range from 88° to 166° F (32° to 74° C)
- Treatment sessions up to 60 minutes
- Specially designed flannel cover for efficient moisture retention
- One-year warranty



Small



Medium



Standard



Shoulder/Neck

## Stainless Steel Carts



### Features

#### Stainless Steel Carts

Versatile, functional and mobile, these sturdy carts go wherever you need them.

### Ordering Information

Item#	Description
Stainless Steel Carts	
A 4002	Model MB – Specially designed to hold the SS, E-2 and E-1 heating units 20" H x 19" W x 17" D (51 cm x 48 cm x 43 cm)
B 4004	Model MB-T – Accommodates stationary heating, ultrasound and electrotherapy units. 30" H x 19" W x 17" D (76 cm x 48 cm x 43 cm)
C 4018	Model MB-TD – Accommodates stationary heating and electrotherapy units. Drawer provides convenient storage for clinic supplies. 30" H x 19" W x 17" D (76 cm x 48 cm x 43 cm)

### Pneumatic Stools with back

77064	Imperial Blue
77065	Burgundy
77066	Dove Gray
77067	Black


### Pneumatic Stools without back

77060	Imperial Blue
77061	Burgundy
77062	Dove Gray
77063	Black

## Pneumatic Therapy Stools



### Available Colors

-  Black
-  Burgundy
-  Dove
-  Imperial Blue

*This color chart is produced within the limitations of lithographic color printing and is intended as a general guide only.*

### Features

#### Pneumatic Therapy Stools

A high-quality, low-priced therapy stool with round, 3" (8 cm) thick foam upholstered swivel seat. Available with or without back in the four colors listed below. A convenient gas spring makes raising and lowering the stool simple 18" - 22" (46 cm - 56 cm). Maximum weight 350 lbs (158 kg).

## Ordering Information

Item#	Description
-------	-------------

Stabilizer™

9296	Stabilizer Unit Comes complete with display gauge, pressure cell and instruction booklet.
------	--

## Stabilizer™



## Features

## Stabilizer™

The Stabilizer Pressure Biofeedback unit gives valuable information to ensure quality in exercise performance and muscle testing. It allows the clinician and patient to determine if the patient is able to selectively isolate and maintain contractions of the cervical or lumbopelvic core stabilization muscles. The measuring range is 0-200 mmHg analog pressure with an accuracy of +/- 3 mmHg pressure.

- Use with exercise to help stabilize joints
- Improve core stabilization exercises
- Enhance home exercise programs



# Sully® Shoulder Support



## Ordering Information

Item#	Size	Chest	Bicep
Sully® Shoulder Support			
660087-001	X-Small	34"-40"	9"-10½"
660087-002	Small	34"-40"	11"-13"
660087-003	Medium	38"-42"	13"-15"
660087-004	Large	42"-46"	14"-17"
660087-005	X-Large	45"-51"	15"-17"

## Features


### Sully® Shoulder Support

On the outside, elastic straps attach with hook to the neoprene at any point, in any direction, and with any amount of force you choose. You can functionally stabilize, assist or restrict movement according to the specific needs of each athlete.


On the inside, the perforated, breathable, rubber-like neoprene grips the skin so the support moves with the athlete. It's as if the stabilization straps are attached directly to the body.

The Sully comes complete with a highly adaptable, bifurcated neoprene strap, and two additional 4" neoprene straps. Worn by professional and amateur athletes in all sports! Measure biceps in a non-flexed position.


## Multiple Applications




Posterior Instabilities




Muscle Strains (pectoral etc.)




Rotator Cuff Deceleration




Anterior Instabilities




Second Strap




First Strap



Second Strap



First Strap



Second Strap

Multi-Directional Instabilities

## Ordering Information

## Item# Description

## A Ultra-Myossage™ Lotion

4262 1 Gallon (3.8 liter) Bottle

## B Myossage™ Lotion

4210 1 Gallon (3.8 liter) Bottle

## C Utility Straps

631374-001 Utility Strap 48"

631325-001 Utility Strap 72"

630483-001 Utility Strap 96"

## D Nylatex Assortment (1 pack of 9 rolls)

1200 Nylatex Assortment

## D Nylatex Individual (3 rolls per package)

1202 \*2.5" x 18" (6 cm x 46 cm)

1264 2.5" x 24" (6 cm x 61 cm)

1204 \*2.5" x 36" (6 cm x 91 cm)

1206 \*2.5" x 48" (6 cm x 122 cm)

1208 \*4" x 18" (10 cm x 46 cm)

1210 \*4" x 36" (10 cm x 91 cm)

1212 \*4" x 48" (10 cm x 122 cm)

1213 4" x 56" (10 cm x 152 cm)

1214 \*6" x 18" (15 cm x 46 cm)

1216 \*6" x 36" (15 cm x 91 cm)

1218 \*6" x 48" (15 cm x 122 cm)

1219 6" x 56" (15 cm x 152 cm)

\* Included in Nylatex Assortment Pack (1200)

## Gels and Lotions



## Features

**A Ultra-Myossage™ Lotion**

Formulated with aloe vera, this ultra-myossage lotion is ideal for massage or as an ultrasound coupling agent. No menthol, hexachlorophene or polysorbates. Water-soluble and greaseless.

**B Myossage™ Lotion**

Soothing, alcohol-free, non-greasy and non-staining lotion.

**C Utility Straps**

A versatile black webbed strap, available in three lengths. Use as a transfer belt, an extension for the EZ-ON traction harness belts, or as shown to assist manual cervical traction.

**D Nylatex Wraps**

- Ideal for holding cold packs in place
- Several sizes available
- Reusable

## Dura-Stick™ Premium Self-Adhesive Electrodes



### Features

#### Dura-Stick™ Premium Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case unless otherwise noted
- Cloth Backing
- Stainless steel mesh
- Choice wire and molded connectors
- Provides current dispersion with low impedance

## Dura-Stick™ Self-Adhesive Electrodes



### Features

#### Dura-Stick™ Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case
- Carbon silver film layer incorporated into a low profile design
- Cloth backing for increased "stretchability" and memory
- Innovative polymer designed to help improved adhesion and durability

Please read complete instructions for use found on electrode packaging prior to use.

### Ordering Information

Item#	Description
DURA-STICK™ Premium Self-Adhesive Electrodes	
42170	Dura-Stick Premium 1.25" Round
42171	Dura-Stick Premium 2" Round
42172	Dura-Stick Premium 1.5" x 2.5" Oval
42173	Dura-Stick Premium 3" x 5" Oval (20/case) (10 packs of 2)
42174	Dura-Stick Premium 2" Square
42175	Dura-Stick Premium 2" x 3.5" Rectangle
42176	Dura-Stick Premium 2" Square Blue Gel
42177	Dura-Stick Premium 1.5" x 3.5" Rectangle Blue Gel
DURA-STICK Self-Adhesive Electrodes	
42160	Dura-Stick 2" x 4" Oval
42178	Dura-Stick 2" Square
42179	Dura-Stick 2" x 3.5" Rectangle


Note: all quantities are 10 packs of 4 each,  
40 per case, unless otherwise stated

## Ordering Information

Item#	Description
Dura-Stick™ Plus Self-Adhesive Electrodes	
42045	Dura-Stick Plus 2.75" x 5" Rectangle
42159	Dura-Stick Plus 1.5" x 2.5" Oval
42181	Dura-Stick Plus 2" x 3.5" Rectangle
42182	Dura-Stick Plus 2" Round
42183	Dura-Stick Plus 2" Square
42160	Dura-Stick Plus 2" x 4" Oval
DURA-STICK Disposable Electrodes	
42180	Dura-Stick 2.25" x 2.5" Disposable 40/Package

Note: all quantities are 10 packs of 4 each,  
40 per case, unless otherwise stated

## Dura-Stick™ Plus Self-Adhesive Electrodes



## Features

## Dura-Stick™ Plus Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case
- Our most economical electrode.
- Permagem™ Hydrogel to help conform to body contours
- Blue foam backing
- Carbon conductive film incorporated into a flexible spun lace backing
- Flexible and sturdy design for durability and ease of use

## Dura-Stick™ Disposable Electrodes



## Features

## Dura-Stick™ Disposable Electrodes

- 40/Package
- Multi-day electrode for extended wear
- May be worn for 2-5 days as guided by a medical professional
- Solid gel center surrounded by adhesive
- Beige mesh backing, looks like a bandage
- Dispose of when finished

Please read complete instructions for use found on electrode packaging prior to use.



## Conductive Rubber Electrodes



### Features

#### Conductive Rubber Electrodes

- Conductive rubber electrodes
- 3" (8 cm) round
- Available in black or red

### Ordering Information

Item#	Description
Conductive Rubber Electrodes	
72852	Black Rubber Carbon 3" (8 cm) Round (2X)
72853	Red Rubber Carbon 3" (8 cm) Round (2X)







**DJO, LLC | A DJO Global Company**

**T** 800.592.7329    **D** 760.727.1280    **F** 800.242.8329  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
[DJOglobal.com/chattanooga](http://DJOglobal.com/chattanooga)