

Sports Medicine Rehabilitation Solutions



◀ Revolution Wireless™
Electrotherapy System



I-Bresis™

Hybrid Iontophoresis System

The Right Medication, the Right Application

- Medications with positive or negative polarity
- Flexible, time-saving treatment options
- Three treatment modes: I-Bresis, Standard, or Patch-Only

The I-Bresis Patch provides fill-and-go efficiency. Built-in batteries plus dose controller activation allows quick and effective application.



Vectra® Neo

- Two and four channel Electrical Stimulation
- Ultrasound
- EMG Biofeedback
- Multimodality versatility with 5 plug-and-play modules
- EMG-Triggered Stimulation
- Laser Options
- Anatomical reference library, including pathologies
- Touchscreen with intuitive graphic user interface



Revolution Wireless™

- Protocol driven programs with evidence-based waveforms
- 14 Pain Management and Muscle Strengthening Programs including programs for: **trigger points, acute and chronic pain, decreasing edema, muscle spasm relaxation, muscle re-education and muscle strengthening**
- 4 channel option for bi-lateral or co-contraction treatment
- Wireless, portable, easy to use



Intellect® Shortwave Diathermy

- Deep tissue heating over a large body area
- Treats both acute and chronic conditions with the convenience of hands-free
- Full color clinic resource library with applicator placement images, anatomical images and pathologies
- Clinical product with fully adjustable parameters

Vectra® Genisys Transport Laser

- Designed with the demands of the Athletic Trainer in mind
- Low-level laser therapy
- Independent control over all laser parameters
- Real-time dosage feedback
- 17 different laser applicators to choose from, including single laser diodes and cluster probes (LED and laser)



Intellect® Transport Ultrasound, Stim, or Combo

- Lightweight design with battery-powered option
- Clinical indications on board along with user-defined memory positions
- Independent controls for each of the two stim channels
- 1 and 3 MHz ultrasound with interchangeable applicator sizes
- Mobile, table top, or cart usage



Primera™ TENS/NMES

- 8 TENS programs, including the HAN TENS waveform, a combination of High and Low Frequencies in a mini sequence
- 6 NMES programs
- Electrode positioning belt included for low back application that can easily be adapted to other areas of the body



M2 Hydrocollator



M4 Hydrocollator



Hotpacs and Terry Covers

Hydrocollator Heating Units, Hotpacs, and Terry Covers

- Hydrocollators – Durable, high-quality stainless steel stationary and mobile heating units of various sizes
- Moist Heat Hotpacs - Seven different sizes to fit various body contours, providing up to 30 minutes of deep, soothing moist heat
- Terry Covers – For high-volume athletic training room use, Hotpac terry covers range in size, shape, and functionality



5 Section



3 Section



2 Section

Montane® Tables

- 5 section, 3 section, 2 section and Bobath platforms
- Electric high-low adjustment from 19" to 40"
- Multilayer padding with select density to provide ideal resistance for the Athletic Trainer yet comfortable for the patient
- Craftsmanship evident in the ultra-quiet actuator and gas springs for lift and powder-coated, laser-cut steel frame with robust, robotic welds
- 360° access with elevation bar



DJO, LLC | A DJO Global Company

T 800.592.7329 **D** 760.727.1280 **F** 800.242.8329
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOGlobal.com/chattanooga

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.