



Moving Rehabilitation Forward®



**Vectra® Neo**  
Simplified, Powerful Solutions...



# Vectra Neo

## Delivering the Results you Demand

Its intelligent design is clever in its features, usability and clinical technology. An exceptional patient and therapist experience.

Every element of the Vectra Neo has been expertly crafted, allowing the clinician to provide unparalleled utilization of the powerful and thoughtful design, empowering them to provide patients with a comprehensive level of rehabilitation.



# You demand versatility and customization

Vectra Neo responds with plug-and-play modules, user-defined menus, and advanced reporting




## Custom made modality

Build the precise combination for your needs. The Vectra Neo has 4 module bays that accommodate the slide-in, plug-in-play module options: Channel 1-2 stim, Channel 1-2 stim/EMG, Channel 3-4 stim, Laser and Ultrasound.

Each unit is assembled and shipped specifically according to your customized clinical needs.

## Customized Menus, Save Customized Treatments

Build your custom programs within the interface, save favorite or patient-specific treatments and customize the system to precisely fit your practice and your patient profiles.



## Customized Patient Reporting

A high speed USB port allows simple and efficient transfer of patient information for therapy documentation and consultation.




# You Demand evidence-based Treatment Options

A convenient, intuitive graphic interface promptly guides you to the treatment protocol

## Clinical Protocol Setup™

Vectra Neo is simple to use thanks to the Clinical Protocol Setup (CPS), which leads you through the functions of the device and each therapy.



With the power of Vectra Neo, you have a comprehensive reference library providing extensive information on therapy, treatments and pathologies – all at your fingertips. By selecting a patient's indications, a list of appropriate protocols will appear. Not only will this aid in facilitating treatments that are tailored to the needs of each patient, it enables you to enhance your daily clinic efficiency. In addition, all of this is managed within the device or can be exported to a USB flash drive.

## Anatomical Library

Vectra Neo offers a rich anatomical library that illustrates an array of pathologies, making it easier for you to communicate with patients about their condition and educate them on further treatment options both in and outside the practice.



# You Demand Form & Function

Smart, innovative and ergonomic design make Neo the high performance choice



## Graphic User Interface

Full-color, 10.4" LCD touchscreen displays commands, data and imagery in sharp, vivid detail, including a gorgeous anatomical library of classic anatomy images and lifelike pathology renderings.

## Lead Wire Management

The innovative lead wire management system ensures access and organization around the full circumference of the unit. Applicators can be inserted into the right or left side for optimal positioning and easy connection.



## Screen Maneuverability

The screen swivels and tilts to a 90-degree range for enhanced patient viewing, privacy and maximized practitioner accessibility.



## Integrated Applicator Cradles

Durable and functional, Neo allows you to house an ultrasound applicator, laser applicator or both - whatever suits your needs. In addition, these applicator cradles are modular and can be detached from the base or switched.

## Cart or Tabletop Configuration

A completely modular setup allows you to choose among the components that best fit your practice; an ergonomic base unit, a cart with 4 module bays, and 5 available treatment options. The unit has an integrated base with discreet, strong handles for carrying. The high-quality cart is stable, height-adjustable, moves with ease, and includes three roomy storage drawers with sturdy pull tracks.

# You Demand Ease of Use

Vectra Neo gives you confidence and control over your outcomes with an intuitive interface

## Precision at Your Finger tips

Vectra Neo seamlessly guides both experienced and new operators through comprehensive menus and treatment parameters with intuitive software and clear, concise instructions. The easy-to-use

touchscreen launches the full spectrum of features all from the home screen, making navigation to treatments, reference sites, patient history quick and simple.



Electrotherapy



Ultrasound



sEMG



Laser



CPS




Custom Protocols



CRPS



Assign Shortcut  
4



Assign Shortcut  
4



Patient Data



Anatomical Library

## Electrode Placement Guide

A comprehensive electrode placement guide with detailed descriptions and imagery.

# You Demand Enhanced Clinical Technology

Vectra Neo's design and detail is the culmination of extensive research and development to bring you a pioneer product in the field of rehabilitation

## **Multipatient**

The Multipatient function gives you the ability to treat more than one\* patient at any one time.

## **Full Selection of Electrotherapy Waveforms**



\*excluding laser

# You Demand Quality & Reliability

Your complete satisfaction is our highest priority, Quality and reliability are the foundation of Vectra Neo and what stands behind the Chattanooga promise everyday.

## Chattanooga

Trust comes from working together. As a leading brand that is dedicated to helping therapists deliver effective rehabilitation solutions, that's exactly what we try to establish. A cornerstone of our development phase is that we work closely with our customers on Voice-of-Customer exercises that ensure product requirements derive from your needs. We track our performance and build on lessons learned, as we constantly strive to deliver you the highest quality.

## Proudly Assembled in the USA

Vectra Neo is proudly designed, assembled and shipped from our DJO Global headquarters in Vista, California.

Learn more about Vectra Neo and prescribe it for your patients; call for a demonstration today and see how a cutting-edge physical agent modality with a simplified operator interface can enhance your ability to achieve optimal outcomes.

Visit [DJOglobal.com/chattanooga/vectraneo](http://DJOglobal.com/chattanooga/vectraneo)



**DJO, LLC | A DJO Global Company**

**T** 800.592.7329    **D** 760.727.1280    **F** 800.936.6569  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
[DJOglobal.com/chattanooga](http://DJOglobal.com/chattanooga)

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.