

3-in-1 Adapter

SENSING

ANALYSIS

DELIVERY

"Designed specifically for use with multiple sizes of medical tubing"



3-IN-1 ADAPTER OVERVIEW

The **3-in-1 Adapter** is designed specifically for use with multiple sizes of commonly used medical tubing. It provides the user with the option to connect 22mm Large Bore, 15mm or 1/4" Cannula tubing to a device or outlet with DISS fittings. A Two-in-One style is available for High Flow applications (22mm and 15mm connection only).

- Up to three connections in a single reusable connector
- Connects Large Bore (22mm), Small Bore (15mm) & 1/4" Cannula to any DISS fitting
- 2-in-1 option for High Flow applications available (22mm & 15mm)
- Made from machined aluminum

PART NUMBERS

3-in-1 (Single)R219P31
3-in-1 (5 Pack)R219P31-001
3-in-1 (50 Pack)R219P31-002
2-in-1 (Single)R219P32
2-in-1 (5 Pack)R219P32-001
2-in-1 (20 Pack)R219P32-002

ML# 295 REV. C