

# 2-in-1 Muffled Adapters

SENSING

ANALYSIS

DELIVERY

*"Dampens noise to minimize patient stimulation during use"*



## PLASTIC 2-IN-1 FEATURES:

- Easily Attach 15mm and 22mm Tubing
- 2-Year Life
- Durable Molded Plastic
- Cleanable
- Low-Cost

## METAL 2-IN-1 FEATURES:

- Easily Attach 15mm and 22mm Tubing
- Extended Long-Life
- Durable Machined Aluminum
- Cleanable



Easily Attach  
15mm and 22mm  
Tubing



## 2-IN-1 MUFFLED ADAPTER OVERVIEW

The **2-in-1 Muffled Adapters** are designed specifically for use with multiple sizes of commonly used medical tubing. Both versions provide users with the option to connect 22mm large bore or 15mm small bore tubing to any device with a male DISS fitting such as common flow meters. They are specially designed for high flow applications and dampen noise to minimize patient stimulation during use.

## PART NUMBERS

### PLASTIC 2-IN-1:

2-in-1, Single . . . . . R219P50-100  
2-in-1, 20 Pack . . . . . R219P50-102  
2-in-1, 50 Pack . . . . . R219P50-101

### METAL 2-IN-1:

2-in-1 Metal, Single. . . . . R219P50

ML# 0248 REV. C