


SMZ160 Series

# SMZ-160-BP





SMZ160 Series

**SMZ-160-BP**

<b>Optical system</b>	Greenough
<b>Observation tube</b>	Binocular head
<b>Inclination</b>	45° inclined
<b>Interpupillary distance</b>	55-75mm
<b>Diopter adjustment</b>	On both tubes, +/- 5 diopter
<b>Eyepieces</b>	Widefield WF10X/20mm
<b>Objectives system</b>	Zoom, ratio 6:1
<b>Objectives</b>	0.75X~4.5X
<b>Working distance</b>	100mm
<b>Stand type</b>	Pole type stand
<b>Head holder</b>	For Ø 25mm column and Ø76mm head
<b>Focus mechanism</b>	Coarse focusing system with tension adjustment
<b>Focusing stroke</b>	43mm
<b>Accessories included</b>	Dust cover
<b>Dimensions LxWxH</b>	280x220x408mm
<b>Base</b>	280x220mm
<b>Column/Arm</b>	Height 215mm, Ø25mm
<b>Net weight</b>	3.2kg
<b>Notes</b>	Transport and storage case